

# ASSURANCE PRODUCTS

BY

ANODYNE SURGICAL


## 0000 SERIES: Retrobulbars and Peribulbars

### Atkinson

0004H Retrobulbar Needle (Atkinson)  
23g x 38mm (1 1/2")

0003H Retrobulbar Needle (Atkinson)  
25g x 38mm (1 1/2")


0001H Retrobulbar Needle (Atkinson)  
27g x 38mm (1 1/2")



0004F Peribulbar Needle (Atkinson)  
23g x 32mm (1 1/4")

0003F Peribulbar Needle (Atkinson)  
25g x 32mm (1 1/4")


0001F Peribulbar Needle (Atkinson)  
27g x 32mm (1 1/4")



0004C Peribulbar Needle (Atkinson)  
23g x 22mm (7/8")


0003C Peribulbar Needle (Atkinson)  
25g x 22mm (7/8")

0001C Peribulbar Needle (Atkinson)  
27g x 22mm (7/8")



0003A Peribulbar Needle (Atkinson)  
25g x 16mm (5/8")

0004A Peribulbar Needle (Atkinson)  
23g x 16mm (5/8")




### Non Atkinson


0008 Retrobulbar Needle (Non Atkinson)  
25g x 38mm (1 1/2")

0009 Retrobulbar Needle (Non Atkinson)  
26g x 38mm (1 1/2")


0010 Retrobulbar Needle (Non Atkinson)  
27g x 38mm (1 1/2")



0011 Peribulbar Needle (Non Atkinson)  
25g x 25mm (1")




0012 Peribulbar Needle (Non Atkinson)  
25g x 22mm (7/8")



**0000 SERIES: Retrobulbars and Peribulbars (cont.)**


**Non Atkinson (cont.)**

- 0013 Peribulbar Needle (Non Atkinson)  
23g x 16mm (5/8")
- 0014 Peribulbar Needle (Non Atkinson)  
25g x 16mm (5/8")
- 0015 Peribulbar Needle (Non Atkinson)  
27g x 16mm (5/8")
- 0016 Peribulbar Needle (Non Atkinson)  
30g x 16mm (5/8")




**Sub Tenon Anesthesia Cannulas**

- 0020 Sub Tenon Anesthesia Cannula  
19g x 25mm (1") curved, flattened tip




- 0021 Sub Tenon Anesthesia Cannula  
21g x 22mm (7/8") curved, tri port

- 0023 Sub Tenon Anesthesia Cannula  
19g x 25mm (1") curved, tri port



- 0022 Sub Tenon Anesthesia Cannula  
21g x 22mm (7/8") 7mm angled tip, tri port

- 0024 Sub Tenon Anesthesia Cannula  
19g x 25mm (1") 7mm angled tip, tri port



- 0027S Rous Sub Tenon/Parabulbar Needle  
25g x 16mm (5/8") Extra Sharp  
curved (bevel down)

The Rous Sub-Tenon/Parabulbar Needle was developed by Dr. Stanley Rous of Fort Lauderdale, FL. It is a 25g x 16mm (5/8") needle gently curved with the bevel on the inside of the curve. This enables safe injection of anesthetic, steroid or antibiotics under direct visualization. The gentle curve and bevel down configuration significantly reduces the risk of potential needle complications compared to long retrobulbar or peribulbar needles that require blind passage and that cannot follow the curvature of the globe.

