











ASSURANCE PRODUCTS


BY

ANODYNE SURGICAL

1000 SERIES: Anterior Chamber Cannulas

	3mm	4mm	5mm	6mm	7mm	8mm	9mm	10mm	11mm	12mm
										
30g	1000-3	*1000-4	1000-5	*1000-6	1000-7	1000-8	1000-9	1000-10	1000-11	1000-12
27g	1001-3	*1001-4	1001-5	*1001-6	1001-7	*1001-8	1001-9	1001-10	*1001-11	1001-12
25g	1003-3	1003-4	1003-5	*1003-6	1003-7	*1003-8	1003-9	1003-10	1003-11	1003-12
23g	1004-3	1004-4	1004-5	*1004-6	1004-7	*1004-8	1004-9	1004-10	1004-11	1004-12
21g	1006-3	1006-4	1006-5	1006-6	1006-7	1006-8	1006-9	1006-10	1006-11	1006-12
20g	1007-3	1007-4	1007-5	1007-6	1007-7	1007-8	1007-9	1007-10	1007-11	1007-12
19g	1008-3	1008-4	1008-5	1008-6	1008-7	1008-8	1008-9	1008-10	1008-11	1008-12

* Denotes most requested Anterior Chamber Cannulas.



Anterior Chamber Cannulas (also known as: Air Injection, Rycroft, and Viscoelastic Cannulas)

Andoyne Surgical offers and stocks the widest variety of Anterior Chamber Cannulas, and has marked the most popular with an "*". To use the chart, select the desired gauge and length of angled tip, then record that catalog number. (For example, a 27g cannula with a 6mm angled tip would be catalog number 1001-6). The drawings at the top of each column depict actual length of the angled tip.

Anterior Chamber Cannulas have the following specifications:

- angle - 32°
- length - 22mm (7/8") PLEASE NOTE: 19g and 20g cannulas are 25mm (1")
- tumble polished tips for smoother radii

NOTE: **Thin Wall** Anterior Chamber Cannulas can be found on pages 6 & 7.

ASSURANCE PRODUCTS

BY

ANODYNE SURGICAL

1000 SERIES: Anterior Chamber Cannulas (Cont.)

1100S Anterior Chamber Cannula (straight)
30g x 22mm (7/8")

1101S Anterior Chamber Cannula (straight)
27g x 22mm (7/8")

1103S Anterior Chamber Cannula (straight)
25g x 22mm (7/8")


1104S Anterior Chamber Cannula (straight)
23g x 22mm (7/8")

1104SL Anterior Chamber Cannula (straight)
23g x 32mm (1 1/4")

1106S Anterior Chamber Cannula (straight)
21g x 22mm (7/8")

1107S Anterior Chamber Cannula (straight)
20g x 25mm (1")

1108S Anterior Chamber Cannula (straight)
19g x 25mm (1")



1100C Anterior Chamber Cannula (curved)
30g x 22mm (7/8")

1101C Anterior Chamber Cannula (curved)
27g x 22mm (7/8")


1103C Anterior Chamber Cannula (curved)
25g x 22mm (7/8")

1104C Anterior Chamber Cannula (curved)
23g x 22mm (7/8")


1106C Anterior Chamber Cannula (curved)
21g x 22mm (7/8")

1107C Anterior Chamber Cannula (curved)
20g x 25mm (1")

1108C Anterior Chamber Cannula (curved)
19g x 25mm (1")



1200 Viscoelastic Cannula (formed)
27g x 22mm (7/8"), formed 11mm from tip



NOTE: Thin Wall Anterior Chamber Cannulas can be found on pages 6 & 7.