



3100 SERIES: Button Coaxial I/A System

- 3100S Button Coaxial I/A System (straight)
17g thin wall outer cannula,
21g inner cannula, 30mm length,
0.3mm side port
- 3100A Button Coaxial I/A System (angled)
17g thin wall outer cannula,
21g inner cannula, 30mm length,
angled 30°,
0.3mm side port




4000 SERIES: "J" Shaped Cannulas


- 4000 "J" Shaped Cannula
30g x 22mm (7/8"),
"J" shape 1.5mm high x 1.5mm wide
- 4001 "J" Shaped Cannula
27g x 22mm (7/8"),
"J" shape 1.5mm high x 1.5mm wide
- 4003 "J" Shaped Cannula
25g x 22mm (7/8"),
"J" shape 1.5mm high x 1.75mm wide
- 4004 "J" Shaped Cannula
23g x 22mm (7/8"),
"J" shape 2mm high x 2mm wide
- 4006 "J" Shaped Cannula
21g x 22mm (7/8"),
"J" shape 2mm high x 2mm wide



- 4003M Micro "J" Shaped Cannula [Hessburg]
25g x 22mm (7/8"),
"J" shape 0.85mm long x 1.7mm wide



- 4003-10 "J" Shaped Cannula
25g x 22mm (7/8"),
"J" shape 2mm high x 2mm wide
10mm angled tip (45°)
- 4004-10 "J" Shaped Cannula
23g x 22mm (7/8"),
"J" shape 2mm high x 2mm wide
10mm angled tip (45°)



ASSURANCE PRODUCTS


BY

ANODYNE SURGICAL

4000 SERIES: "J" Shaped Cannulas (Cont.)

4003-17L "J" Shaped Cannula (angled left) [Binkhorst]
25g x 22mm (7/8"),
"J" shape 1.5mm high x 1.75mm wide
12mm angled tip (18°)

4003-17R "J" Shaped Cannula (angled right) [Binkhorst]
25g x 22mm (7/8"),
"J" shape 1.5mm high x 1.75mm wide
12mm angled tip (18°)



5000 SERIES: Cystotomes

Irrigating Cystotomes

5000AS Irrigating Cystotome (straight)
30g x 16mm (5/8")

5001AS Irrigating Cystotome (straight)
27g x 16mm (5/8")


5001BS Irrigating Cystotome (straight)
27g x 19mm (3/4")

5003AS Irrigating Cystotome (straight)
25g x 16mm (5/8")

5003CS Irrigating Cystotome (straight)
25g x 22mm (7/8")

5004AS Irrigating Cystotome (straight)
23g x 16mm (5/8")

5004BS Irrigating Cystotome (straight)
23g x 19mm (3/4")



5000AF Irrigating Cystotome (formed)
30g x 16mm (5/8")

5001AF Irrigating Cystotome (formed)
27g x 16mm (5/8")


5001BF Irrigating Cystotome (formed)
27g x 19mm (3/4")

5003AF Irrigating Cystotome (formed)
25g x 16mm (5/8")


5003CF Irrigating Cystotome (formed)
25g x 22mm (7/8")

5004AF Irrigating Cystotome (formed)
23g x 16mm (5/8")


5004BF Irrigating Cystotome (formed)
23g x 19mm (3/4")



16mm (5/8")



19 mm (3/4")



22mm (7/8")