

ASSURANCE PRODUCTS


BY

ANODYNE SURGICAL

4000 SERIES: "J" Shaped Cannulas (Cont.)

4003-17L "J" Shaped Cannula (angled left) [Binkhorst]
25g x 22mm (7/8"),
"J" shape 1.5mm high x 1.75mm wide
12mm angled tip (18°)

4003-17R "J" Shaped Cannula (angled right) [Binkhorst]
25g x 22mm (7/8"),
"J" shape 1.5mm high x 1.75mm wide
12mm angled tip (18°)



5000 SERIES: Cystotomes

Irrigating Cystotomes

5000AS Irrigating Cystotome (straight)
30g x 16mm (5/8")

5001AS Irrigating Cystotome (straight)
27g x 16mm (5/8")


5001BS Irrigating Cystotome (straight)
27g x 19mm (3/4")

5003AS Irrigating Cystotome (straight)
25g x 16mm (5/8")

5003CS Irrigating Cystotome (straight)
25g x 22mm (7/8")

5004AS Irrigating Cystotome (straight)
23g x 16mm (5/8")

5004BS Irrigating Cystotome (straight)
23g x 19mm (3/4")



5000AF Irrigating Cystotome (formed)
30g x 16mm (5/8")

5001AF Irrigating Cystotome (formed)
27g x 16mm (5/8")




5001BF Irrigating Cystotome (formed)
27g x 19mm (3/4")

5003AF Irrigating Cystotome (formed)
25g x 16mm (5/8")

5003CF Irrigating Cystotome (formed)
25g x 22mm (7/8")

5004AF Irrigating Cystotome (formed)
23g x 16mm (5/8")

5004BF Irrigating Cystotome (formed)
23g x 19mm (3/4")



ASSURANCE PRODUCTS


BY

ANODYNE SURGICAL

5000 SERIES: Cystotomes (Cont.)


5001AA Irrigating Cystotome (angled)
27g x 16mm (5/8")

5003AA Irrigating Cystotome (angled)
25g x 16mm (5/8")




5001AC Irrigating Cystotome [Joyce]
27g x 16mm (5/8") curved

5003AC Irrigating Cystotome [Joyce]
25g x 16mm (5/8") curved




5003ASC Irrigating Cystotome (side cut)
25g x 16mm (5/8")




5001ASR Irrigating Cystotome (small radius)
27g x 16mm (5/8")

5003ASR Irrigating Cystotome (small radius)
25g x 16mm (5/8")



5003ARF Irrigating Cystotome (reverse formed)
25g x 16mm (5/8")
tip angled 90° towards bevel




Capsulorhexis Cystotomes

5007 Capsulorhexis Cystotome (formed)
27g x 16mm (5/8")
tip twisted 90°

5007S Capsulorhexis Cystotome (straight)
27g x 16mm (5/8")
tip twisted 90°

5008 Capsulorhexis Cystotome (formed)
25g x 16mm (5/8")
tip twisted 90°



ASSURANCE PRODUCTS

BY


ANODYNE SURGICAL

5000 SERIES: Cystotomes (Cont.)

5008S Capsulorhexis Cystotome (straight)
25g x 16mm (5/8")
tip twisted 90°

5009 Capsulorhexis Cystotome (formed)
23g x 16mm (5/8")
tip twisted 90°


5009S Capsulorhexis Cystotome (straight)
23g x 16mm (5/8")
tip twisted 90°



Irrigating Sharp Cystotomes

5010 Irrigating Sharp Cystotome (straight)
[Kelman] 22g x 25mm (1")


5011 Irrigating Double Sharp Cystotome (straight)
[Kelman] 22g x 25mm (1")



6000 SERIES: Capsule Polishers

6000 Capsule Polisher (side port)
21g x 22mm (7/8")
7mm angled tip (37°)
3mm sandblasted closed distal end
0.3mm side port

6001 Capsule Polisher (side port)
23g x 22mm (7/8")
7mm angled tip (37°)
3mm sandblasted closed distal end
0.3mm side port



6002 Capsule Polisher (open end)
21g x 22mm (7/8")
7mm angled tip (37°)
3mm sandblasted distal end

6003 Capsule Polisher (open end)
23g x 22mm (7/8")
7mm angled tip (37°)
3mm sandblasted distal end

