



8000 SERIES: Miscellaneous Instruments


8000 Lacrimal Cannula (curved)
25g x 28mm (1 1/8")
tumble polished edge




8000S Lacrimal Cannula (straight)
19g tapered to 23g 12mm from end



8001 Lacrimal Intubation Set
stainless steel probes
23g x 175mm (7")
silicone tubing length
300mm (12")
tubing dimensions
0.3mm (.012") x 0.6mm (.025")



8002 Lacrimal Intubation Set
stainless steel probes
23g x 110mm (4 1/4") with bulbous tips
silicone tubing length
300mm (12")
tubing dimensions
0.3mm (.012") x 0.6mm (.025")
retriever is stainless steel with a 0.75mm opening
on a ribbed autoclavable plastic handle.
handle length 75mm (3")




* Ribbed non-slip handle

8002WO Lacrimal Intubation Set
(no retriever)

8003S DCR Set (straight)
stainless steel probes 20g x 45mm (1 3/4")
silicone tubing length 300mm (12")
tubing dimensions 0.5mm (.020") x 0.95mm (.037")


8003A DCR Set (angled)
stainless steel probes 20g x 45mm (1 3/4")
silicone tubing length 300mm (12")
tubing dimensions 0.5mm (.020") x 0.95mm (.037")



8000 SERIES: Miscellaneous Instruments (Cont.)


8100 Stiletto Knife 20g
1.15mm width
handle length 75mm (3")

8101 Stiletto Knife 19g
1.38mm width
handle length 75mm (3")




8200S Micro Lens Hook (straight) [Sinskey II]
used for 0.2mm holes or larger
handle length 75mm (3")


8200A Mirco Lens Hook (angled) [Sinskey II]
10mm angled tip (45°)
used for 0.2mm holes or larger
handle length 75mm (3")




8201 Lens Manipulating Hook [Fenzl]
10mm angled tip (45°)
used for 0.2mm holes or larger
handle length 75mm (3")




8202 Button Tip Manipulator (angled)
10mm angled tip (45°)
handle length 75mm (3")



8203 Lens Manipulator [Jaffe - Maltzman] (angled)
10mm angled tip (45°)
"V" shaped blunt tip
handle length 75mm (3")



8204 IOL Manipulating Hook [Lester] (angled)
10mm angled tip (45°)
hourglass shaped tip
handle length 75mm (3")




All instruments on this page have ribbed non-slip handles

8000 SERIES: Miscellaneous Instruments (Cont.)


8205 Nucleus Rotator [Jaffe-Bechert] (angled)
7mm angled tip (45°)
"Y" shaped blunt tip
handle length 75mm (3")

* Ribbed non-slip handle




8206 Micro Iris Retractor [Pearce] (angled)
1mm wide tip

* Ribbed non-slip handle




8207 Irrigating Dialer
30g x 22mm (7/8")
10mm angled tip (45°)
1mm distal end angled 65°




8208 Nucleus Spatula [Pearce]
.80mm wide X .20mm thick
10mm angled tip (45°)

* Ribbed non-slip handle




8209 Nucleus Chopper [P. Koch] (angled)
6mm angled tip (45°)
handle length 75mm (3")

* Ribbed non-slip handle




8400 Eye Speculum (adult)[Barraquer]

8401 Eye Speculum (child) [Barraquer]



8402 Eye Speculum (adult) [Kratz]
open blades



ASSURANCE PRODUCTS

BY

ANODYNE SURGICAL

8000 SERIES: Miscellaneous Instruments (Cont.)


8500 Alpha-chymotrypsin Cannula (curved)
25g x 28mm (1 1/8")
1mm smooth, bulbous tip

8500S Olive Tip Capsule Polisher (curved)
25g x 28mm (1 1/8")
1mm sandblasted bulbous tip


8501 Olive Tip Cannula (curved)
23g x 22mm (7/8")
1mm smooth, bulbous tip

8501S Olive Tip Capsule Polisher (curved)
23g x 22mm (7/8")
1mm sandblasted, bulbous tip


8503 Olive Tip Capsule Polisher
23g x 10mm bend (angled)
1mm sandblasted, bulbous tip




8600 Infusion Handle
75mm (3"), luer slip end and tapered tubing end
made with autoclavable plastic




8601 Male/Female Hand Piece
75mm (3"), fits luer taper on both ends
made with autoclavable plastic



8602 Double Male Hand Piece
50mm (2"), fits luer taper on both ends
made with autoclavable plastic



8603 Aspirating Hand Piece with finger hole
75mm (3"), fits luer taper on both ends
made with autoclavable plastic



9000 SERIES: Chamber Maintainers and Retinal Cannulas

9000 Chamber Maintainer
27g x 22mm (7/8")
0.2mm side port, sharp closed end
6" silicone tubing: ID - 0.3mm (0.012")
OD - 0.6mm (0.025")

9001 Chamber Maintainer
25g x 22mm (7/8")
0.2mm side port, sharp closed end
6" silicone tubing: ID - 0.3mm (0.012")
OD - 0.6mm (0.025")

