



anodynesurgical.com  
 804 Corporate Centre Drive  
 O'Fallon, Missouri, USA 63368

Phone: (800) 428-5628

*Vision in Motion*

## II DISPOSABLE KNIVES

### SLIT KNIVES



#### Single Bevel Slit Knife

IMAGE	PRODUCT NO	DESCRIPTION
	AE02-1801-4506	1.8mm Single Bevel Angled Slit Knife
	AE02-2001-4506	2.0mm Single Bevel Angled Slit Knife
	AE02-2201-4506	2.2mm Single Bevel Angled Slit Knife
	AE02-2401-4506	2.4mm Single Bevel Angled Slit Knife
	AE02-2501-4506	2.5mm Single Bevel Angled Slit Knife
	AE02-2601-4506	2.65mm Single Bevel Angled Slit Knife
	AE02-2701-4506	2.75mm Single Bevel Angled Slit Knife
	AE02-2801-4506	2.8mm Single Bevel Angled Slit Knife
	AE02-2901-4506	2.85mm Single Bevel Angled Slit Knife
	AE02-3001-4506	3.0mm Single Bevel Angled Slit Knife
	AE02-3201-4506	3.2mm Single Bevel Angled Slit Knife

#### Double Bevel Slit Knife

IMAGE	PRODUCT NO	DESCRIPTION
	AE02-1803-4506	1.8mm Double Bevel Angled Slit Knife
	AE02-2203-4506	2.2mm Double Bevel Angled Slit Knife
	AE02-2403-4506	2.4mm Double Bevel Angled Slit Knife
	AE02-2503-4506	2.5mm Double Bevel Angled Slit Knife
	AE02-2603-4506	2.6mm Double Bevel Angled Slit Knife
	AE02-2703-4506	2.75mm Double Bevel Angled Slit Knife
	AE02-2803-4506	2.8mm Double Bevel Angled Slit Knife
	AE02-2903-4506	2.85mm Double Bevel Angled Slit Knife
	AE02-3003-4506	3.0mm Double Bevel Angled Slit Knife
	AE02-3203-4506	3.2mm Double Bevel Angled Slit Knife

## II DISPOSABLE KNIVES

### SLIT KNIVES




#### Clear Cornea Slit Knives

IMAGE	PRODUCT NO	DESCRIPTION
	AE02-1804-4506	1.8mm Angled Clear Cornea Slit Knife
	AE02-2204-4506	2.2mm Angled Clear Cornea Slit Knife
	AE02-2404-4506	2.4mm Angled Clear Cornea Slit Knife
	AE02-2504-4506	2.5mm Angled Clear Cornea Slit Knife
	AE02-2604-4506	2.6mm Angled Clear Cornea Slit Knife
	AE02-2704-4506	2.75mm Angled Clear Cornea Slit Knife
	AE02-2804-4506	2.8mm Angled Clear Cornea Slit Knife
	AE02-3004-4506	3.0mm Angled Clear Cornea Slit Knife
	AE02-3204-4506	3.2mm Angled Clear Cornea Slit Knife

### STAB KNIVES



IMAGE	PRODUCT NO	DESCRIPTION
	AE01-1503-0006	15 deg Stab Knife
	AE01-2203-0006	22.5 deg Stab Knife
	AE01-3003-0006	30 deg Stab Knife
	AE01-4503-0006	45 deg Stab Knife

[www.anodynesurgical.com](http://www.anodynesurgical.com)

804 Corporate Centre Drive O'Fallon, Missouri, USA 63368

(800) 428-5628



Anodyne products are used in 46 countries and growing



Anodyne products are used in 5 million (+) cases per year



Privately-owned company established in 1981

## II DISPOSABLE KNIVES

### MVR STYLE KNIVES



#### Sideport Knives

IMAGE	PRODUCT NO	DESCRIPTION
	AE01-1003-4506	1.0mm Angled MVR Style
	AE01-1103-4506	1.1mm Angled MVR Style
	AE01-1203-4506	1.2mm Angled MVR Style

#### MVR Blades

IMAGE	PRODUCT NO	DESCRIPTION
	AE04-1903-4506	19g Angled MVR Knife
	AE04-1903-0006	19g Straight MVR Knife
	AE04-2003-4506	20g Angled MVR Knife
	AE04-2003-0006	20g Straight MVR Knife
	AE04-2303-4506	23g Angled MVR Knife
	AE04-2303-0006	23g Straight MVR Knife
	AE04-2503-4506	25g Angled MVR Knife
	AE04-2503-0006	25g Straight MVR Knife

### CRESCENT KNIVES



#### Crescent Knives

IMAGE	PRODUCT NO	DESCRIPTION
	AE03-2601-4506	2.6mm Angled Crescent Bevel Up
	AE03-2001-4506	2.0mm Angled Crescent Bevel Up
	AE03-2601-0006	2.6mm Straight Crescent Bevel Up

[www.anodynesurgical.com](http://www.anodynesurgical.com)

804 Corporate Centre Drive O'Fallon, Missouri, USA 63368

(800) 428-5628